
SPIRO®

HUMIDIFIER

REF 702 CE 0402

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sv Bruksanvisning

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Symbol Glossary

Symbol	Title	Description
	Caution!	Indicates the need for the user to consult the Instructions for Use for important cautionary information, such as warnings and precautions.
	Consult instructions for use	Indicates the need for the user to consult the Instructions for Use.
	Manufacturer	Indicates the medical device manufacturer.
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Use-by date	Indicates the date after which the medical device is not to be used.
	Batch code	Indicates the manufacturer's batch code so the batch or lot can be identified
	Keep away from sunlight	Indicates a medical device that needs protection from light sources.
	Keep dry	Indicates a medical device that needs to be protected from moisture.

Symbol	Title	Description
	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	Do not re-use	Indicates a medical device that is intended for one use. (Article 740 only)
	Single patient only	Indicates that the device is only to be used by a single patient.
	Medical device	Indicates that the device is a medical device.
	Quantity	Indicates the total number of devices in the package. X is a variable.
Rx Only	Prescription only (USA)	Indicates that the product is a medical device as defined in 21 CFR 820.3(1) and Federal Law (USA) restricts this device to sale by or on the order of a physician (21 CFR 801.109)
	CE mark	Indicates the manufacturer declaration that the product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.
	Distributor	Indicates the distributor of the devices and their contact details for re-ordering.

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https://ec.europa.eu/health/md_sector/contact_el

SPIRO HUMIDIFIER (REF 702)Single patient use only 

INTENDED USE: The Humidifier is intended for use by spontaneous breathing tracheostomy patients. When placed on the tracheal cannula the Humidifier uses the exhaled air to warm and moisten the inhalation air. The Humidifier has a removeable COVER (5) to allow the HME (heat moisture exchange) filter to be changed. The Humidifier is intended for both adult and paediatric populations. The Humidifier has an a 15-mm-ISO-connector and shall be replaced with a new unit after 29 days of use, or more frequently if necessary.

N.B. THE HUMIDIFIER MUST NEVER BE USED WITHOUT THE FILTER

DAILY CHECKS: Check humidifier operation daily and inspect for damage. Replace humidifier with new unit in case of damage or defective operation.

CLEAN the valve daily. Do not use sharp tools. Use the fingers or **PLASTIC TWEEZERS** (REF 705, sold separately) to remove/re-insert the filter. The valve **HOUSING** (3) and **COVER** (5) shall be air dried after washing with luke-warm water and mild fragrance-free detergent. Other detergents or solvents must not be used. Always check the function again after cleaning and check that the air passages are free. **DO NOT WASH AND RE-USE THE FILTER!**

REPLACE THE FILTER (4) (REF 740, sold separately) at least twice every 24 hours, or more frequently if required. Replacement may be required more often, perhaps several times per hour, if a build-up of mucus or phlegm is present. Ensure that the filter is seated evenly against the bottom of the filter housing, and that the air flow through the filter is not obstructed by phlegm or other substances.

A **SHOWER PROTECTOR/HEAT SHIELD** (6) (REF 703, sold separately) should be used as required with the inhalation opening facing downward.

An **ADAPTOR** (2) (REF 707, 708, sold separately) is only required for connection to a cannula (1) without a 15-mm-ISO connector (taper).

The valve's accessories (SHOWER PROTECTOR, ADAPTORS and TWEEZERS) are reusable after daily cleaning according to the valve's cleaning instructions above, and shall be replaced with a new unit after 29 days of use.

The Humidifier **REF** 702 is supplied in BASIC SET SPEAKING VALVE **REF** 760/760R, and BASIC SET HME **REF** 770. Each set is supplied with a Speaking Valve (**REF** 760/760R only), Shower Protector for Humidifier **REF** 703, Tweezers **REF** 705, Adaptors **REF** 707/708, and Shower Protector for Speaking Valve **REF** 730 (**REF** 760/760R only).

SUPPLEMENTAL OXYGEN can connect via the Ø5.5mm nipple on the valve housing (Max. 5 litres/minute).

RESISTANCE TO FLOW for SPIRO Humidifier is < 0.8 cm H₂O and with heat shield/shower protection on it has a resistance to flow < 1.2 cm H₂O (at 30 litres per minute).

The Humidifier must be removed before inhalation of medicine.

All products should be stored in a dark and dry environment and at normal indoor temperature.

For the STORAGE PERIOD of SPIRO PRODUCTS: See year and month next to the 'HOURGLASS' symbol on the product label. The products shall be disposed of in accordance with the applicable national regulations.

ALWAYS CONSULT YOUR DOCTOR, SUPPLIER OR MANUFACTURER IF YOU HAVE ANY QUESTIONS.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

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SPIRO Basic Device Description

REF	Description	Basic UDI-DI
702	SPIRO Humidifier	7350061230016V
703	SPIRO Shower Protector for Humidifier (REF 702)	7350061230026X
705	SPIRO Tweezers	7350061230036Z
707	SPIRO Adaptor 13/11 mm	73500612300473
708	SPIRO Adaptor 14/11 mm	73500612300473
710	SPIRO Speaking Valve with oxygen nipple (transparent)	73500612300575
710R	SPIRO Speaking Valve with oxygen nipple (red)	73500612300575
711	SPIRO Speaking Valve (transparent)	73500612300677
730	SPIRO Shower Protector for Speaking Valve (REF 710/710R/711)	73500612300779
740	SPIRO Replacement Filter for Humidifier (REF 702) and Speaking Valve (REF 710/710R/711)	7350061230087B
760	SPIRO Basic Set Speaking Valve; containing Speaking Valve (REF 710), Humidifier (REF 702), Shower Protector for both Humidifier (REF 703) and Speaking Valve (REF 730), Tweezers (REF 705), Adaptors (REF 707+708)	7350061230097D

SPIRO Basic Device Description

REF	Description	Basic UDI-DI
760R	SPIRO Basic Set Speaking Valve; containing Speaking Valve ([REF]710R), Humidifier ([REF]702), Shower Protector for both Humidifier ([REF]703) and Speaking Valve ([REF]730), Tweezers ([REF]705), Adaptors ([REF]707+708)	7350061230097D
770	SPIRO Basic Set HME; containing two Humidifiers ([REF]702), Shower Protector for Humidifier ([REF]703) Tweezers ([REF]705), Adaptors ([REF]707+708)	7350061230106W

Fogless International AB
Storgatan 24
SE-567 32 Vaggeryd
Sweden
Tel. +46-36-144100
Fax. +46-36-144116
mail@fogless.se
www.fogless.se