

SPIRO[®]

SPEAKING VALVE

REF 710 CE 0402

REF 710R CE 0402

REF 711 CE

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sv Bruksanvisning

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Symbol Glossary

Symbol	Title	Description
	Caution!	Indicates the need for the user to consult the Instructions for Use for important cautionary information, such as warnings and precautions.
	Consult instructions for use	Indicates the need for the user to consult the Instructions for Use.
	Manufacturer	Indicates the medical device manufacturer.
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Use-by date	Indicates the date after which the medical device is not to be used.
	Batch code	Indicates the manufacturer's batch code so the batch or lot can be identified
	Keep away from sunlight	Indicates a medical device that needs protection from light sources.
	Keep dry	Indicates a medical device that needs to be protected from moisture.

Symbol	Title	Description
	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	Do not re-use	Indicates a medical device that is intended for one use. (Article 740 only)
	Single patient only	Indicates that the device is only to be used by a single patient.
	Medical device	Indicates that the device is a medical device.
	Quantity	Indicates the total number of devices in the package. X is a variable.
Rx Only	Prescription only (USA)	Indicates that the product is a medical device as defined in 21 CFR 820.3(1) and Federal Law (USA) restricts this device to sale by or on the order of a physician (21 CFR 801.109)
	CE mark	Indicates the manufacturer declaration that the product complies with the essential requirements of the relevant European health, safety and environmental protection legislation.
	Distributor	Indicates the distributor of the devices and their contact details for re-ordering.

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https://ec.europa.eu/health/md_sector/contact_el

SPIRO SPEAKING VALVE **REF** 710 + 710R + 711Single patient use only 

INTENDED USE: The Speaking Valve is a one-way valve intended for use by spontaneous breathing tracheostomy patients. THE SPEAKING VALVE ALLOWS FOR INHALATION BUT NOT EXHALATION THROUGH A TRACHEAL CANNULA. The valve opens upon inhalation and closes upon exhalation, channelling the expired air up through the vocal cords and mouth/nose. The Speaking Valve is intended for both adult and paediatric populations.

WARNING: NEVER PLACE A SPEAKING VALVE ON A TRACHEAL CANNULA THAT HAS AN INFLATED CUFF AND IS NOT FENESTRATED. THERE MUST BE AN OPEN PASSAGE FOR THE EXHALATION AIR TO PASS UP THROUGH THE VOCAL CORDS WHEN USING A SPEAKING VALVE.

CONTRAINDICATIONS: This device is not for use with patients who have respiratory infection, extensive phlegm, stenosis, vocal cord paralysis or acute airway obstructions. This device is not for use with inline ventilators, endotracheal tubes or other artificial airways.

PRIOR TO USE: The Speaking Valve connects to a tracheal cannula with a 15-mm-ISO-connector. Ensure that the patient is awake and alert. Ensure that the patient has a tracheal cannula with a deflated cuff or no cuff; and/or fenestration. Make sure the patient can exhale through the mouth/nose before placing the Speaking Valve on the patient's cannula. Measurement of the tracheal pressure at passive exhalation with the Speaking Valve in place is advised (< 5 cm H₂O is recommended).

The Speaking Valve has a removeable COVER (5) to allow for changing of the filter. The cover has a tapered section which seats in the valve HOUSING (3) during use. The further pressed in, the tighter the cover is secured. Avoid pressing the cover too hard, as it reduces the valve's life and makes it difficult to remove.

DAILY CHECKS: Check the function of the valve daily and inspect for damage. Replace valve with new unit in the case of damage or defective operation. **NEVER DISASSEMBLE THE MEMBRANE!**

CLEAN the valve daily. Do not use sharp tools. Use the fingers or PLASTIC TWEEZERS (**REF** 705, sold separately) to remove/re-insert the filter. The valve housing and cover shall be air dried after washing with luke-warm water and mild fragrance-free detergent. Rinse the valve cover in both directions to remove all cleaning liquids. Other detergents or solvents must not be used. Always check the operation again after cleaning and check that the air passages are free. DO NOT WASH AND RE-USE THE FILTER!

The valve should always be used with FILTER (4) (**REF** 740, sold separately). REPLACE THE FILTER daily or more frequently if required. Ensure that the filter is seated evenly against the bottom of the valve housing, and that air flow through the filter is not obstructed by phlegm or other substances.

A SHOWER PROTECTOR/HEAT SHIELD (6) (**REF** 730, sold separately) should be used as required with the inhalation opening facing downward.

An ADAPTOR (2) (**REF** 707, 708, sold separately) is only required for connection to a cannula (1) without a 15-mm-ISO-connector (taper).

The valve's accessories (SHOWER PROTECTOR, ADAPTERS and TWEEZERS) are reusable after daily cleaning according to the valve's cleaning instructions above, and should be replaced with a new unit after 29 days of use.

Speaking Valve **REF** 710/710R has a Ø5.5 mm nipple for SUPPLEMENTAL OXYGEN on the valve housing (3) (Max. 5 litres/minute). The WHITE O2 NIPPLE CAP (7) must always be fitted to the Speaking Valve and only removed when oxygen is being supplied.

Speaking Valve **REF** 710/710R is supplied in BASIC SET SPEAKING VALVE **REF** 760/760R respectively. Each set is supplied with Humidifier **REF** 702, Shower Protector for Humidifier **REF** 703, Tweezers **REF** 705, Adaptors **REF** 707/708, and Shower Protector for Speaking Valve **REF** 730.

The SPEAKING VALVE shall be replaced with a new unit after 29 days of use, or more frequently if required. The valve must not be used while the patient is sleeping or otherwise when the patient's consciousness is impaired. Use a SPIRO HUMIDIFIER (**REF** 702, sold separately) in such cases.

The Speaking Valve must be removed before inhalation of medicine.

RESISTANCE TO FLOW for SPIRO Speaking Valve is < 2 cm H₂O, and with the heat shield/shower protection on < 3 cm H₂O (at 30 litres per minute).

All products should be stored in a dark and dry environment and at normal indoor temperature.

For the STORAGE PERIOD OF SPIRO PRODUCTS: See year and month next to the 'HOURGLASS' symbol on the product label. The products shall be disposed of in accordance with the applicable national regulations.

ALWAYS CONSULT YOUR DOCTOR, SUPPLIER OR MANUFACTURER IF YOU HAVE ANY QUESTIONS.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

The contacts of national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

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SPIRO Basic Device Description

REF	Description	Basic UDI-DI
702	SPIRO Humidifier	7350061230016V
703	SPIRO Shower Protector for Humidifier (REF 702)	7350061230026X
705	SPIRO Tweezers	7350061230036Z
707	SPIRO Adaptor 13/11 mm	73500612300473
708	SPIRO Adaptor 14/11 mm	73500612300473
710	SPIRO Speaking Valve with oxygen nipple (transparent)	73500612300575
710R	SPIRO Speaking Valve with oxygen nipple (red)	73500612300575
711	SPIRO Speaking Valve (transparent)	73500612300677
730	SPIRO Shower Protector for Speaking Valve (REF 710/710R/711)	73500612300779
740	SPIRO Replacement Filter for Humidifier (REF 702) and Speaking Valve (REF 710/710R/711)	7350061230087B
760	SPIRO Basic Set Speaking Valve; containing Speaking Valve (REF 710), Humidifier (REF 702), Shower Protector for both Humidifier (REF 703) and Speaking Valve (REF 730), Tweezers (REF 705), Adaptors (REF 707+708)	7350061230097D

SPIRO Basic Device Description

REF	Description	Basic UDI-DI
760R	SPIRO Basic Set Speaking Valve; containing Speaking Valve ([REF] 710R), Humidifier ([REF] 702), Shower Protector for both Humidifier ([REF] 703) and Speaking Valve ([REF] 730), Tweezers ([REF] 705), Adaptors ([REF] 707+708)	7350061230097D
770	SPIRO Basic Set HME; containing two Humidifiers ([REF] 702), Shower Protector for Humidifier ([REF] 703) Tweezers ([REF] 705), Adaptors ([REF] 707+708)	7350061230106W

Fogless International AB
Storgatan 24
SE-567 32 Vaggeryd
Sweden
Tel. +46-36-144100
Fax. +46-36-144116
mail@fogless.se
www.fogless.se